
TENSIOMETRE POIGNET EXTRA PLAT

FR

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REF : T2021 - MD2210

CONSIGNES DE SECURITE

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
DIGITAL AUTOMATIC WRIST BLOOD PRESSURE MONITOR

USER MANUAL

REF : T 2021 - MD2210

SAFETY NOTES

WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

WARNING

- Consult your physician before starting to measure your blood pressure.
- Like any oscillometric blood pressure measurement devices, certain medical conditions can affect the measurement accuracy, among others:
 - Disorder of the cardiac rhythm
 - Very low blood pressure
 - Very low blood perfusion
 - Patients in shock
 - Diabetes
 - Vessel anomalies
 - People with electrical implants such as a cardiac pacemaker
 - Women who are pregnant
 - Patient who are pre-eclamptic
 - Patient motion during measurement

Due to their condition the oscillometric measurement method can produce incorrect readings. This represents a risk for your health, since values may be interpreted incorrectly. Always consult your physician to determine what will be suitable for you.

SAFETY NOTES

- This product does not and is not intended to provide a medical diagnosis. Measurement results are for reference only. Self-diagnosis and treatment, e.g. regarding medication, using measured results represent a risk for your health. Always consult with a licensed physician for determination of appropriate medication and dosage thereof. Follow the instructions of your physician or licensed healthcare provider.
- Please note that technically related measuring tolerances are possible. Please see the section “TECHNICAL SPECIFICATION”.
- Any cuff related blood pressure measurement in high repetition rates, can lead to severe measurement side effects, e.g.
 - A nerve compression with temporary wrist/ hand paralysis
 - The release of an arterial or venous thrombus, which can cause a life threatening situation.Please contact your physician about the specific risks of cuff pressure in your specific case.
- The “irregular heartbeat” function does not replace a cardiac examination, but may help to detect potential heart rate irregularities at an early stage. Always consult your physician to determine what will be suitable for you.
- The “irregular heartbeat” function is not designed for diagnosing or treating an arrhythmic disorder. Arrhythmia can only be ascertained by a licensed physician.
- The “WHO BLOOD PRESSURE CLASSIFICATION” chart is not intended to replace a medical diagnosis. This chart is only for reference for different classifications of blood pressure.
- If you notice abnormal or suspicious variations in blood pressure measurements, consult your physician immediately.
- Women who underwent a breast or axillaries lymph node removal operation should consult a physician or licensed healthcare provider before starting blood pressure measurements.
- This device must be used in accordance to the specified ambient conditions, otherwise the accuracy of readings might be affected. See “TECHNICAL SPECIFICATION” for details.

SAFETY NOTES

- Do not wrap the cuff around body parts other than your left wrist. Misuse represents a risk to your health.
- This product is not intended for use by or on children, toddlers and infants or on persons who cannot express their consent, e.g. persons with mental disorders or the like. Consult your physician for alternative methods of measuring a child's blood pressure.
- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- This product contains small parts that may present a choking hazard to children. Keep the unit and all parts out of reach of children.
- Proper cuff size is important for accurate measurements. Only use the device on adults who have the right wrist circumference for this unit. See "TECHNICAL SPECIFICATION" for suitable wrist circumferences.
- Electromagnetic interference: Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile phones, microwave ovens) while it is in operation, as inaccurate measurements may result. To prevent such interference, use the unit at a sufficient distance from such devices or turn the disturbances off.
- Batteries should not be charged or reactivated by any other means. The batteries may explode.
- Take extra precaution to keep a leaking battery away from fire as there is a risk of ignition or explosion.
- In case battery fluid leaks and come into contact with your eyes or skin, do not rub and immediately rinse with plenty of clean water and seek medical advice.
- Do not use the equipment where flammable gas (e.g. anaesthetics gas, oxygen or hydrogen) or flammable liquid (e.g. alcohol) are present.

SAFETY NOTES


CAUTION

- This device should not be used when your wrist has been wounded/ injured or when a catheter has been inserted. Such use may result in injury.
- Remove any kind of wrist jewellery or the like before taking a measurement. This could cause bruises.
- Do not place the wrist cuff over heavy clothing (e.g. a jacket or sweater sleeve) as the blood pressure monitor will not be able to take a proper measurement and there is an elevated danger of acquiring hematoma or skin marks during the course of the measurement.
- In case the cuff does not stop inflating, interrupt the measurement by pressing the ON/OFF button and open the cuff at once.
- In case battery fluid leaks, do not touch the battery fluid. Avoid skin contact (e.g. put on protective gloves) and clean the battery compartment with dry cloth.
- Do not disassemble the device, it may cause damage to yourself and to the product. If you cannot fix the problem using the "TROUBLESHOOTING" instructions, request service from your dealer.
- When applying the cuff, make sure there are no wrinkles in the cuff as this could cause bruises.
- Blood pressure measurements can lead to temporary marks on the skin at the site of the cuff placement. This is especially the case in high repetition rates, in hypertonic patients and in patients with weak heart rates, In rare cases a mark may persist for couple of days. Please contact your physician about these specific risks of cuff pressure in your specific case.

MAINTENANCE AND STORAGE


- The blood pressure monitor is constructed by precision electronic components, Accuracy of readings and the instrument's service life depend on careful handling. Protect the unit against hard knocks (e.g. dropping the unit), moisture, water, dirt, dust, chemicals, extreme hot or cold temperatures, major temperature fluctuations, direct exposure to sunlight and heat sources which are too close (e.g. stoves, heating radiators). This may damage the unit. The device must be stored in the specified ambient conditions. Please see section "Technical Specification" for details.
- The device is designed and manufactured for a long service life. However it is generally recommended to have the monitor inspected every 2 years to ensure proper functioning and accuracy. Please contact your dealer for maintenance.
- Never immerse and/or spill water or any other liquid onto the monitor or any components, otherwise liquid may enter it and cause damage.
- Never use rechargeable batteries. This may damage the unit.
- Replace all batteries at the same time and use batteries of the same type. Do not mix old and new batteries.
- Never attempt to repair, open and/or disassemble the unit or adjust it yourself. This may cause damage to the unit and impair functions. If you cannot fix the problem using the "TROUBLESHOOTING" instructions, request service from your dealer.
- Do not drop or insert any object into any openings or hoses. This may damage the unit.
- Do not press the buttons with excessive force or with pointed objects.
- Clean your device and cuff carefully only with a slightly moistened soft cloth and dry it immediately with a soft dry cloth. Do not press.
- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals (e.g. thinner, alcohol, benzene) to clean the device.
- When storing the device, make sure that no heavy objects are placed on top of it.
- Do not fold the cuff and tubing tightly. The cuff tube should not have any sharp kinks and keep it away from sharp edges.
- Leaking batteries may damage the unit. If you do not intend to use the unit for longer periods, remove the batteries from the battery compartment before placing the device in storage.
- Used equipment, parts and batteries are not treated as ordinary household waste, and must be disposed according to applicable local regulations for material disposal. Unlawful disposal may cause environmental pollution.

PART IDENTIFICATION & DISPLAY READINGS


INSTALL/ CHANGE BATTERIES


1. Pull the battery cover located on the bottom side on the unit as shown in figure.

2. Remove old batteries and insert new ones. Use LR03 or AAA alkaline batteries and use batteries of same type only.
3. Make sure the battery polarities (+) and (-) match the marking on the battery compartment.

4. Place back the battery door.

NOTE:

- Insert the batteries as shown in the battery compartment. If not, the device will not work or even be damaged.
- [] icon represent battery power is enough.
- When battery power becomes weak, the [] icon and "E6" appears in the display, replace all batteries with new ones. Do not mix old and new batteries.
- Never use rechargeable batteries. This may damage the unit.
- Battery life may vary with ambient temperature and may be shorter at low temperature.


TO SET CLOCK

When new batteries are installed

1. "YEAR" will blink on display automatically.
2. Press [+ / M] to select current Year.
3. Press [⏻] button to confirm and "MONTH" will display.
4. Press [+ / M] to select current Month.
5. Press [⏻] button to confirm and "DAY" will display.
6. Press [+ / M] to select current Day.
7. Press [⏻] button to confirm and "HOUR" will display.
8. Press [+ / M] to select current Hour.
9. Press [⏻] button to confirm and "MINUTE" will display.
10. Press [+ / M] to select current Minute.
11. Press [⏻] button to confirm and clock setting is done.
12. The set date and time will be displayed.

At any time

1. [⏻] button to when the device is off.
2. Date and time will be displayed.
3. Press and hold [⏻] button for 3 seconds and "YEAR" will blink.
4. Follow the same procedure above to set clock.


TO APPLY WRIST CUFF

1. Bare your left wrist, make sure that the circulation in the arm is not restricted by any cloth etc. that are too tight. Position the cuff on the inside of your wrist.

2. Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approximate 1 cm below the ball of your thumb.

3. The cuff has to be fitted tightly around the wrist but should not constrict it.

BEFORE MEASURING BLOOD PRESSURE


1. Apply the wrist cuff following the instruction in “TO APPLY WRIST CUFF”.
2. Rest at least for 5 minutes before each measurement. Otherwise there may be erroneous results.
3. Sit down in a comfortable position, which your legs should be uncrossed, feet should be feet on the floor and your back should be supported by a chair. At the same time your wrist should be supported on a flat surface which the cuff is on a level with your heart.
4. Relax your wrist and turn your palm upward.
5. Relax, keep still and do not talk during the measurement in order not to distort the result.

NOTE:

- This instrument is intended for use by adult only. Do not use this device on or by children, toddlers and infants.
- For reliable monitoring and reference of blood pressure, it is recommended to do the measurement daily at the same time.
- In order to obtain a resting condition blood pressure, do not eat, drink alcohol and caffeinated beverages, smoke, exercise or bath for at least 30 minutes before taking a measurement since your blood pressure varies time to time depending on what you have eaten, drunk and what you have done.
- To minimize measurement variations due to physical activity, relax for five to ten minutes before taking a measurement. If you are excited by emotional stress or muscular movement, the measurement result given out may be inaccurate. Stress raises blood pressure.
- You should not be physically tired or exhausted while taking measurement.
- Perform measurements in a quiet and relaxed environment at room temperature.
- Remain relax, still and do not speak during measurement since the accuracy of any blood pressure measurement with this device can be affected by those multitude of causes.
- Always wait at least 5 minutes between measurements to allow the blood circulation in your wrist to return to normal. You may need to increase the wait time depending on your individual physiological characteristic.
- Should the device detect an abnormal condition, it stops the measurement and display an error code. See “TROUBLESHOOTING” for more details.
- This device bases it measurements on your heartbeat. If you have a very weak or an irregular heartbeat, the device may have difficulty in determining your blood pressure. (An irregular heartbeat is defined as a heartbeat that varies by 25% from the average of all heartbeats during measurement.)

TO MEASURE YOUR BLOOD PRESSURE

1. Follow the instruction in "BEFORE MEASURING BLOOD PRESSURE" to make yourself ready.
2. Press [①] button when the device is off to start measurement. The average blood pressures and last data number are displayed briefly. The current pressure is displayed, as the measurement begins.
3. The cuff starts to inflate. It is normal for the cuff to feel very tight. A pressure bar indicator is displayed during measurement. See "PRESSURE BAR INDICATOR" for more details.
4. When inflation is complete, deflation starts automatically. Once the pulse is detected, the [♥] blinks with each pulse beat, indicating that the measurement is in progress.
5. When the measurement is complete, the systolic and diastolic pressure and pulse rate are displayed and stored. At the same time, the cuff exhausts the remaining air and deflates completely.
6. Press [①] button to turn off the power.


Note:
You can stop inflation by pressing [①] button at any time.

PRESSURE BAR INDICATOR

Standards to assess high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart below.


Reference Material: *Journal of Hypertension* 1999, Vol 17 No.2


WHO CLASSIFICATION INDICATOR

Each of the six segments of the bar indicator corresponds to the WHO blood pressure classification.

WHO Classification Indicator:

Example:

The indicator displays a segment, based on the current data, corresponding to the WHO classification.

TO RECALL AVERAGE AND PREVIOUS MEASUREMENT DATA

This device has a memory capability to store the measurement readings for each user. Every time you complete the measurement, the device automatically stores the measurement result.


1. To get the average measurement data of the last 3 measurements, press [] key to turn on the unit.
2. Press [] key to select the user you desire.
3. Press [+/M] to read the average data.
4. Press [+/M] continuously to view previous measurement data.
5. Press [] button to turn of the device.

TO DELETE MEASUREMENT DATA

This monitor has 2 methods to delete the recorded memory.


(A) To delete one single record

1. Select the record you want to delete
2. Press and hold [+/M] button until the screen display "CL" and "--"
3. Press [] button to confirm and the screen will display "CL" and "00"

(B) To delete all memories

1. Press the [+/M] button to enter average mode
2. Press and hold [+/M] button until the screen display "CL" and "--"
3. Press [] button to confirm and the screen will display "CL" and "00"

WHAT IS AN IRREGULAR HEARTBEAT

This blood pressure monitor provides a blood pressure and pulse rate measurement even when an irregular heartbeat occurs. An irregular heartbeat is defined as a heartbeat that varies by 25% from the average of all heartbeats during the blood pressure measurement. It is important that you are relaxed, remain still and do not talk during measurements.

Note: We recommend contacting your physician if you see this [♥👉] indicator frequently.

ABOUT BLOOD PRESSURE

What Is Blood Pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands. Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

What Is Hypertension And How Is It Controlled?

Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering one's lifestyle, avoiding stress, and with medication under a doctor's supervision. To prevent hypertension or to keep it under control:

- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

TROUBLESHOOTING

Nothing appears in the display, even when the power is turned on	Batteries are drained	Replace all batteries with new ones
	Battery polarities are not in the correct position	Re-install the batteries with their negative and positive ends matching their indicated in the battery compartment
ERROR code 1 (E1) appears	The cuff position is not fastened properly	Fasten the cuff correctly
	The cuff position is not correct	Sit comfortably and still. Ensure that the cuff is the same level as the heart
ERROR code 2 (E2) appears	You moved your arm or body during measurement	Make sure you remain very still and quiet during the measurement

TROUBLESHOOTING

ERROR code 3 (E3) appears	The cuff may not be applied	Check whether tube connection of the cuff is secured to the unit properly
ERROR code 4 (E4) appears	The unit does not measure	If you have a very weak or irregular heart beat, the device may have difficulty in determining your blood pressure
	There is a measuring Error	Sit comfortably and still. Fasten the cuff again carefully
ERROR code 5 (E5) appears	Cuff over inflated	The measurement range is over 300 mmHg. It is recommended to see doctor as soon as possible.
ERROR code 6 (E6) appears	Low battery	The battery power is too low to function. Replace the batteries with new ones.
The monitor keeps reinflating	Circuit locked	Remove and reinsert the batteries and then proceed to take measurement again.

TECHNICAL SPECIFICATION

Display	: LCD Display
Measurement Range	: Systolic Pressure: 50-250 mmHg ; Diastolic Pressure: 30-200 mmHg Pulse : 40-180 beats/minute
Accuracy	: Pressure : +/-3 mmHg Pulse : +/-5% of reading
Resolution	: Pressure : 1 mmHg Pulse : 1 beat / minute
Measurement Method	: Non-invasive, Oscillometric method
Power Source	: 2 x 1.5V AAA batteries
Operating Temperature / Humidity	: +5°C to +40°C, 15-93% RH maximum
Storage Temperature / Humidity	: -25°C to +70°C, up to 93% RH maximum
Operation, storage and transport atmospheric pressure	: 700hPa to 1060hPa
Outer Dimensions	: Approx. 90 x 68 x 18mm
Wrist Circumference	: 13.5 - 23cm
Accessories	: Cuff, User manual, Storage box, Batteries(Optional)
Classification	: Application part Type BF
Key to symbols	: Application part Type BF  Class II equipment symbol 

SYMBOLS

Symbols	Function / Meaning
SN	Serial Number
	Manufacturer
	Type BF: Device, cuff and tubing are designed to provide special protection against electrical shocks.
SYS	Systolic Blood Pressure in mmHg
DIA	Diastolic Blood Pressure in mmHg
PUL	Pulse
 0123	EC Directive Medical Device Label
	Caution
	Authorized Representative in the European Community
	WEEE Label
	Refer to instruction manual / booklet
	Keep dry

NOTICE:

- Operation mode : Continuous

The European directive 2012/19/EC about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capacities, or by persons lacking in experience or knowledge, unless they are supervised or instructed on the use of this apparatus by a person responsible for their safety. Children should be supervised to ensure they do not play with this appliance.

Precision and precautions about the battery use

- Throwing away a used battery in the environment or with household waste pollutes and avoids the recuperation of reusable materials.
- It is important to limit the battery consumption and to respect the following rules:
 - Favor alkaline battery (last longer than saline battery) and when it's possible use rechargeable battery.
 - Leave used batteries and accumulators in adapted containers that you can find in some retailers. In that way, metals don't pollute the environment and the health.
- Battery has to be inserted by respecting the polarity show on the apparel and the battery. An incorrect positioning can damage the apparel, cause leaking and even provoke fire and the explosion of the battery.
- To assure a good functioning, battery has to be in a good condition. If there is a problem, put new battery.
- Never try to recharge battery not made for this. It could leak, warm up, cause fire or explode.
- Replace all the battery at the same time. Never mix up alkaline battery with saline battery.
- Used battery has to be removed.
- Also, remove the battery if you're not using the apparel for a long time.
- Never try to short-circuit the terminal of a battery.
- Never throw away the battery in fire because they can explode.
- Take out the accumulator from the apparel before reload it.

- *We recommend that an adult watch if a child change the battery to respect those rules.*
- *If you swallow a battery, go immediately to your doctor or the hospital.*
- *For the best ongoing performance, we recommend using only fresh, high-quality, alkaline batteries.*

Importé par / Importiert von / Geïmporteerd door / Imported by :
CAPTELEC, 59170 Croix - FRANCE

Carton et emballage papier à
trier / Kartons und
papierverpackung recyceln /
Verpakkingsmateriaal in karton
en papier te sorteren / Recycle
carton and paper packaging