94247

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## GB - BMI tape measurer

The BMI is a measurement of your body weight based on your height and weight. Although your BMI does not actually "measure" your percentage of body fat, it is a useful tool to estimate a healthy body weight based on your height. Due to its ease of measurement and calculation, it is the most widely used diagnostic indicator to identify a person's optimal weight depending on his height. Your BMI "number" will inform you if you are underweight, of normal weight, overweight, or obese. However, due to the wide variety of body types, the distribution of muscle and bone mass, etc., it is not appropriate to use this as the only or final indication for diagnosis.

### Directions for use:

The small central rotary circle indicates the size. The great outdoors and fixed circle indicates the weight. Turn the mobile circle to place your height in front of your weight and your BMI you can read on the scale window:

Between 10 and 20: You are underweight

Between 20 and 25: You are normal

25 to 30: You are overweight

Between 30 and 40: You are obese

Between 40 and 60: You are very obese

You can also take measurements of different parts of your body through the meter built into the device. Press the button on the back of the machine to wind meter.