DE

Akupressurmatte (ca. 73 x 42 x 2 cm, 500g)

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EN

Acupressure mat (approx. 73 x 42 x 2 cm, 500 g)

Bring a piece of Asian medicine and wellness culture into your home.

Acupressure is an ancient East Asian healing method, in which a therapeutic effect shall be achieved by means of thin needles that are inserted into specific points of the body. In traditional Chinese medicine (TCM), the application of acupuncture is based on "life energy of the body" (Qi).

The perfect all-rounder for the treatment of tension and blockages as well as for the stimulation of your blood circulation. A total of 8820 pressure tips, distributed on 210 "needle flowers" (made of ABS), comfortably massage your back, abdomen and feet. After you have been using it, you will be surrounded by a pleasant tranquility and pain-relieving relaxation.

Instructions for use:

Lay down gently on the mat. Depending on how much experience you already have with the usage, mild pain symptoms may initially occur after three to five minutes. Listen to your body's signals when using it. If you feel too much unpleasant, put a piece of cloth between the mat and your body and dress yourself lightly. Within the first 14 days, you can increase the dosage from 10 to 20 minutes. The evening hours are particularly well suited due to the fact that a pleasant and slightly tedious feeling enhances the relaxation character.

Different positions allow the acupressure to be applied to different parts of the body.

Position 1: The back

Place the mat on a flat surface and lay down flat on it. You can also put a pillow under your neck in order that you can lay more comfortable. In this position, the acupressure points stimulate along the spine in order to prevent any back tensions.

Position 2: The stomach

Lay on your stomach and cross your arms under your head. This position helps to improve metabolism and also supports the relaxation of breathing and diaphragm.

Position 3: The jaw

Place your cheek on the rolled-up acupressure mat. This position relaxes jaw and neck muscles.

Depending on your wishes and needs, you can also adopt these 3 positions one after another or mixed.

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