
**Pedaltrainer RFM<sup>®</sup> digital**

**DEUTSCH**

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
Pedaltrainer Exerciser RFM® digital

ENGLISH


- 1) Please take the **Pedal Exerciser RFM® digital** carefully out of its' transport box (see picture 1)
- 2) Unfold the stands into their locked position, place the **Pedal Exerciser RFM® digital** on a smooth surface and screw in the adjustment lever A (see picture 2)
- 3) Adjust the pressure on axis B by turning the lever A. Turning the lever left increases the pressure, while turning it right decreases the pressure. Turn it back to Minimum after use. Press button C to activate the built-in digital unit, repeated pressing switches trough the programs (see picture 3)
- 4) Pull the safety pins to fold and stock the Pedal Exerciser RFM digital after use (see pictures 4 and 5)
- 5) Please note:
  - (a) Check the unit for damages before using it. Make sure the stands are locked and the exerciser is places on firm ground.
  - (b) Readjust the pressure for axis B before use as per your requirements.

## The Digital Unit:

The displays:

<p><b>Time</b></p> <p>Press C 1x:                    arrows at CNT and Time</p> <p>The unit shows total revolutions and time, which continues to tick on even when you stop pedalling!</p>	<p><b>Caloric Expenditure</b></p> <p>Press C 2x:                    arrows at CNT and CAL</p> <p>The unit depicts total revolutions and caloric expenditure.</p>
<p><b>RPM</b></p> <p>Press C 3x:                    arrows at CNT and RPM</p> <p>The units shows total revolutions and RPM</p>	<p><b>Scan Mode</b></p> <p>Press C 4x:                    arrows at CNT, RPM and Scan</p> <p>The unit shows total revolutions and time of use, caloric expenditure and RPM (last 3 values in intervals of about 10 seconds).</p>

Unit will switch off after about 3 minutes without training

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