Artikel
Item
Article

Artikel
Articolo
Articulo

4681

4681GAz1_1308

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1 | Introduction

Thank you for choosing this warming product. This is a quality product that should last a long time when used properly.

The "Rechargeable Warming Belt" is supplied with power by two Lithium Ion rechargeable batteries which are integrated in the device's controller. The degree of warmth can be selected in 4 levels between 38 and 55 degrees Centigrade.

2 | Important safety hints

Retain the manual for further use!

- Forward the manual to other users.
- This device is not intended to be used unsupervised by children, individuals with sensory deprivation of any kind or mentally challenged individuals. Any users falling into these categories should be assisted by a responsible adult.

Caution!

- Do not prick in with needles or other sharp objects.
- Do not fold blanket if still warm.
- Do not place heavy or edged objects on the blanket, which may damage or bend the electrical wires inside the blanket (eg. Suit cases).
- Do not place other heaters or heat producing objects onto the switched on device.

3 | Intended use

This product is designed to be an auxiliary source of warmth for individual use that is unencumbered by wires tied to a power source. Sensation of the warming effect is dependent on the clothing of the user. Therefore wearing thick coats or other heavy clothing, the sensation of the warmth will be delayed.

The technical safety of this product is assured by its compliance to european

standards. It is however not intended to treat medical conditions

4 | Start-up

Please charge the battery fully before using the device for the first time. To accomplish this, plug the charger **5** into a power mains and connect with the controller/battery **2** and allow charging for 3-4 hours. The extent of the charge is indicated by a red LED **7** on the controller. The indicator light will turn green once charge is complete.

Afterwards: Remove charger. **5**

4.1 | Operating instructions

- a Connect the cable of the warming belt **1** with the controller. **2**
- b To turn on: Depress the Start/Stop key **6** for two seconds.
- c Select the desired level of warmth by pressing the Start/Stop key **6** repeatedly. It has 4 levels of heat.
- d Place the controller **2** in the side pocket **3** during use.
- e To turn off: depress and hold the Start/Stop key **6** for two seconds. The LED **7** which indicated the level of warmth will turn off.

Caution!

A red blinking LED light, shows a charge is required.

5 | Cleaning and storing

- Use a soft, slightly damp cloth to wipe.
- The Massage Belt is not machine-washable.
- When cleaning, do not submerge in liquids.
- Store in a clean, dry place.

6 | Delivery Contents

- Rechargeable warming cushion
- Controlling electronics(controller) and Lithium Ion rechargeable battery
- Mains adaptor
- Belt extender

- Pocket for controller with hook and loop closure
- Instruction Manual

Not included!

The battery/controller component of this device can be procured under the part number 4681.1.33 in case a second battery is desired. This controller also matches to our rechargeable warming belt art 4681.

7 | Technical Data

Model/Type:	LH-115
Mains adaptor:	100-240 V~ 50-60 Hz / 8,4 V
Input voltage:	8,4 V
Power:	9 W
Material:	100% Polyamid
Weight:	110 g
Product dimensions:	ca. 100 x 16,5 cm

8 | Disposal

Electrical and electronic devices may not be disposed of with domestic waste. Consumers are obliged by law to return electrical and electronic devices at the end of their service lives to the public collecting points set up for this purpose or point of sale. Details to this are defined by the national law of the respective country. This symbol on the product, the instruction manual or the package indicates that a product is subject to these regulations. By recycling, reusing the materials or other forms of utilising old devices, you are making an important contribution to protecting our environment.

9 | Warranty

The guarantee lasts for 24 months from the date of purchase for manufacturers' error in material and quality. If it is used for another purpose or operated incorrectly, no liability for possible damages can be accepted. This guarantee is only valid if the instructions in the "instructions for use" have been followed.

10 | Returns

Dear Customer,

In case you product does not function as you expected please do not reurn immediately. Perhaps our Hotline (+49-69-954061-24) can be of assistance If you are still not satisfied, you have the possibility to return products under 40 € value free of charge within the BRD. We require your address. This you may provide to the above listed telephonic Hotline or by email to : technik@hydass.de.

Prepare the package for shipping and include the receipt and description of the problem. Please be aware of the fact that collect packages may not be accepted.

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