WiTouch Pro

Quick Start Guide

Congratulations on your purchase of the Hollywog® WiTouch Pro device. Hollywog WiTouch Pro is the first wireless remote controlled pain relief device incorporating TENS technology to specifically target back pain. The thin and flexible design perfectly contours the back for maximum surface contact. The advanced electronics design maximizes energy use, providing over 150 30-minute treatment sessions per battery life. This innovative device is safe, drug-free, easy to use, discreet & comfortable to wear, and most importantly allows you to control your pain to maintain an active lifestyle.

Get started in a few easy steps

This Quick Start Guide gives you the basics to use your WiTouch Pro. It will introduce you to the features and explain how to get started with your pain relief treatment.

Learn more

Please read the User Manual for complete information on your WiTouch Pro. The User Manual includes detailed information on intended use, contraindications, precautions, warnings, operation and maintenance of your WiTouch Pro. For additional support, visit www.hollywog.com.

Caution: Federal law [USA] restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device. Use and application of this device is to be considered and advised by the licensed practitioner prior to prescribing the device.

Get to know your WiTouch Pro

What comes with your WiTouch® Pro?

- 1 Hollywog® WiTouch Pro device
- 1 Hollywog® WiTouch Pro remote control
- 1 Package Hollywog® WiTouch gel pads (5 pair)
- 2 AAA alkaline batteries
- 1 CR2032 lithium battery (installed in remote)
- 1 Size #0 Phillips head screwdriver tool
- 1 Storage Box
- 1 User Manual
- 1 Quick Start Guide

Operate your WiTouch Pro

1. Prepare the treatment area by washing and drying the skin thoroughly. Excessive hair should be trimmed (NOT shaved), and the treatment area should be void of oils and/or lotions. DO NOT apply the main device over broken skin.
2. Remove battery cover from WiTouch Pro device with included screwdriver tool and insert (2) AAA batteries noting correct +/- polarity and replace cover.

3. Press and hold the ON/OFF button on the WiTouch Pro device for one (1) second. The LED will begin blinking green once every two (2) seconds. Your WiTouch Pro device is now ON and ready for use.
4. Remove clear battery insulation tab from the remote control to activate battery. Discard clear tab.

5. Remove two (2) WiTouch gel pads from package and separate the pads at the perforation in the liner.

6. Remove blue liners from gel pads by peeling away the folded edges of the liners. Align the shape of the gel pad to match the shape of the black electrode area, and apply gel pads to WiTouch Pro device. Firmly press the gel onto the electrode across the entire surface of the gel pad to ensure good adhesion.

DO NOT remove green liners during this step.

7. Remove green liners from gel pads by slowly peeling the liner from one inside corner towards the outside corner of the electrode area and apply WiTouch Pro device on the back over area of pain, aligning the center of the unit with the spine.

DO NOT activate stimulation of the device prior to application of the device to the back.

8. Press and hold START/STOP button on remote control momentarily to start treatment stimulation mode. On the WiTouch Pro device, the LED will now begin blinking green very rapidly for the duration of the treatment.

9. Press INCREASE INTENSITY (+) button to increase stimulation intensity to a comfortable level.

NOTE: The stimulation treatment will automatically stop after 30 minutes. To start a new treatment, repeat step 8. The Hollywog WiTouch Pro device will turn off after two (2) hours of inactivity. To turn on the device, repeat step 3.

Hollywog

2830 Amnicola Hwy
Chattanooga, TN 37406 USA
+1 423 305 7778

info@hollywog.com • www.hollywog.com
® and ©2013 Hollywog, LLC. All Rights Reserved
Patents Pending

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© and ©2013 Hollywog, LLC. All Rights Reserved