91241

Imported by - Importé par:  
GD import (Genius Ideas & Wellys)

190 chemin de la Frayère - 06530 Peymeinade - France

www.gdimport.com

www.genius-ideas.com

www.wellys-care.com

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FR- Hydropropulseur "Push and Clean"/ 209025 version 2016-07-12

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GB- POWER FLOSS "Push and Clean"/ 209025

Diagram 1: While using the power floss, do not put it in horizontal position, it might have some leaks.

- A. RESERVOIR
- B. MAX FILL LINE
- C. TUBE
- D. BUTTON

- E. PUMP HOUSING
- F. NOZZLE

#### PREPARATION:

**Important:** Before using the item, disconnect the reservoir and submerge the pump housing, tube and nozzle completely in water.

Press the button for at least 60 seconds until it starts bubbling.

Then lift the end of the nozzle above the water and press the button to ensure that the unit is working properly and has a steady jet stream.

#### FILLING YOUR POWER FLOSS:

1. Hold the power floss upright, gently twist the reservoir counter clockwise and separate from the pump housing.
2. Fill the reservoir up to the max fill line with either water or mouthwash. **NOTE: use mouth wash with the same consistency as water.**
3. Reattach the pump housing by inserting it into the reservoir and gently twist clockwise to lock it.

#### USING YOUR POWER FLOSS:

1. Position your thumb as shown in Diagram 1 on the button.
2. Always keep the item in an upright position during use to ensure a steady jet stream.
3. Press the button several times until liquid comes out in a steady stream. You can adjust the speed and volume of the stream by varying the angle of the power floss as you move it across your mouth.
4. To clean in-between teeth: position the nozzle against the teeth at the gum line, close your mouth around the nozzle and press the button. Use at both, front and back of the teeth for optimum results and continue for as long as required.
5. To clean Bridges, braces, implants or other devices: position the nozzle under it and press the button. For best results, irrigate each site from more than one angle whenever possible.

#### Caring for your power floss:

1. After each use, empty the reservoir and depress the button to remove any remaining fluid from within the unit.
2. For hygienic reasons do not share it.
3. Store it in a safe place and ensure the nozzle cannot be damaged.
4. When travelling, store the nozzle in the reservoir.

#### CLEANING YOUR POWER FLOSS:

After use, empty the reservoir and depress the button to remove any remaining fluid from within the unit.

Flush it thoroughly using warm water and wipe dry.

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