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Z 09841 M DS V5 0818

Dear Customer,

We are delighted that you have chosen this **acupuncture pen**. This article is intended for the symptomatic treatment of acute and chronic musculoskeletal pain (excluding pain resulting from bone fractures). It is also indicated for the symptomatic treatment of migraine. **In the case of medical indication, consult your doctor before use!**

If you have any questions, contact the customer service department via our website: www.service-shopping.de

WARNING – Do not use the article:

- if you use a pacemaker.
- if you suffer from heart rhythm problems or epilepsy.
- if you have previously suffered from skin irritation, itching or reddening of the skin or have experienced any other adverse reaction from using TENS.
- during the first 16 weeks of pregnancy and around the belly area throughout the entire pregnancy. Pregnant women should only use the article following consultation with their doctor.
- if taking anticoagulant or steroid treatments, or in case of special clinical conditions, without first consulting a doctor.
- while driving.
- in the presence of gas, flammable liquids or surroundings rich in oxygen.
- on skin where you have applied creams, oils or any other flammable or alcoholic products.
- on infants or people who are unconscious.

Failure to follow these warnings could result in serious injury.

Precautions:

- Read the instructions and keep them dry and in a safe place.
- Keep out of reach of children.
- The article has not been tested on children under 10 years of age. **The article is not suitable for children under the age of 10!**
- Do not apply the article at the carotid sinus node, on the laryngeal and pharyngeal muscles, near the heart, near the eyes, on the mucous membranes.
- Do not apply to injured or irritated skin.
- **Do not apply directly over metal implants.**
- Do not apply on wet or unclean skin or in damp surroundings. **Do not use the article when it is wet or damp.**
- Do not use if the article is damaged and internal parts are exposed.
- Always hold the article in the correct manner. Firmly grip around the metal ring encircling the article and ensure the tip is in contact with the area to be treated.
- No product or advice is intended as a substitute for a prescription or medical diagnosis. The article is only suitable for external use and not for therapeutic purposes. It only relieves pain and does not replace medical treatment.
- Please consult your doctor if you are concerned about the applicability of the article in relation to your condition.
- Stop using the article immediately if you suddenly feel weak or dizzy, if you suffer from skin irritation/reddening, itching or other side effects.
- **WARNING – Danger of Suffocation!** Keep children and animals away from plastic bags and films.
- Protect the article from: impacts, heat, direct sunlight, frost, moisture and wet conditions.
- No do modify the article. For repairs contact the customer service department.

Explanation of Symbols

Danger symbol: Read through the associated safety instructions and follow them.

Read the instructions prior to use!

Symbol for devices classified as "Type BF" (level of protection against electric shock)

Manufacturer within the meaning of the Medical Devices Act

CE 1282

CE conformity pursuant to the European Directive 93/42 CEE of the medical devices under the supervision of Notified Body number 1282

Technical Data

Article number:	09841
Classification:	Medical device belonging to class IIa (CE 1282)
Energy emitted:	2.6 mJ at 2000 Ω
Type of protection:	IPX1
No voltage supply required.	
Transport and storage:	
	– Temperature: -10°C to +50°C
	– Relative humidity: 10% to 95%

CE 1282

Disposal

The packaging material should be recycled. Dispose of the article in the household waste.

Use on Pain Points

1. Take hold of the acupuncture pen and place your index finger around the metal ring.
2. Place the tip of the pen on the pain point or acupuncture point and press the trigger up to 40 times.

- The pen will also work through light clothing.
- In order to treat larger areas of pain, move the pen over these areas and repeatedly press the trigger as you do so.

Each click produces an electrical pulse that feels like a light pinch. This is intended to stimulate the release of endorphins (the body's own pain-relieving agents).

The electrical pulses are generated by a piezoelectric crystal. This is sufficient for at least 100,000 clicks.

Use on Acupuncture Points

This article can be used for a non-invasive acupuncture treatment. The following pictures show 26 acupuncture points. They are followed by a list of pain areas with numbers for those acupuncture points which are to be treated for the respective pain.

If more than 3 points are listed for one pain area, try the first 3 points first. If the effect is unsatisfactory, choose different points until you have found the best combination. Click on each acupuncture point that you want to treat 20-25 times.

If you are treating a migraine or headache, you must perform the stimulation on the opposite side of the body. If the headache or migraine occurs in the right half of the head, perform the stimulation on the left side of your body, for example on your left arm or on your left hand.

Acupuncture Points

Pain Areas

Head

- Head injury or headache: 8, 3, 7, 6, 2, 10, 1, 11, 20, 23
- Headache:
 - General: 1, 3, 7, 4, 6, 2, 10, 11, 19, 20
 - Migraine: 1, 7, 12
- Neck: 10, 3, 14, 1, 7, 12
- Toothache: 4, 1, 18, 21
- Ears: 3, 1

Trunk

- (Nighttime) leg cramps: 4
- Chest: 3, 2, 1, 7
- Abdomen:
 - Upper abdomen: 5, 6
 - Lower abdomen: 9, 5, 6, 2, 19
- Ovaries and womb: 8, 13, 5, 6
- Hernia: 5, 6, 19
- Sciatica: 9, 4, 5, 6, 24, 22
- Muscle cramps: 25, 26
- Lung: 7
- Stomach: 8, 5, 6, 2
- Menstruation: 13, 5, 2, 15, 26, 19
- Neck: 8, 1, 7, 16, 17, 26, 23
- Back pain:
 - General: 9, 4, 5, 6, 17, 20, 22, 23
 - Shoulder blade region: 17, 22, 23, 9, 4
 - Small of the back: 4, 5, 6, 20
- Muscle pain caused by straining or overstretching a muscle: 4
- Lumbago: 9, 4, 5, 6, 14, 1, 7, 24
- Testicles: 4, 5, 6
- Urinary tract/genitals: 13, 25, 26, 19
- Whiplash: 8, 4, 1, 7
- Neuralgia: 4

Arms and legs

- Arms: 8, 3, 2, 12, 17, 23
- Arthritis/osteoarthritis, gout: Use the acupuncture pen on the pain point
- Legs: 9, 4, 5, 6
- Thighs: 9
- Elbows: 8, 3, 2, 14, 1, 23
- Ankles: 4, 5, 6, 1, 25, 18, 19, 20
- Hands: 3, 2, 10, 1, 12
- Hips: 4, 5, 6
- Knee: 9
- Armpits: 2, 12
- Forearms: 8, 2, 10, 1, 12
- Wrists: 7, 3, 12
- Shoulders: 8, 3, 2, 14, 1, 7, 17, 26
- Fingers: 10, 3, 2, 1, 12
- Feet: 4, 5, 6, 1, 25, 18, 19, 20

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